

PROTOCOL INFORMATION

Disorder: Diabetic Foot Wound Healing Date: May 4, 2011
Developed By: Dr. Hao Wu Reviewed By: Dr. Paul Ziemer

PROTOCOL DESCRIPTION

Disorder Description: Diabetic foot wound healing can be accelerated by increasing temperature, blood supply and local oxygen ratio in operated area. These three factors are critical for tissue recovery process. Genesis light can improve nerve sensitivity in diabetic neuropathy.

Treatment Position: 1. High light 20 min or Low light 30 min over diabetic foot wound



Treatment Protocol: Time: 20mins Position: 1 or more Light: High Treatment Interval: 1-2 time/day, 8-12 weeks

Result: Satisfactory Acceptable Desirable Exceptional Fully healed